
DETAIL LEGEND:

- 1) HUNTER DECODER (ICD-200) PER PLAN
- 2 HUNTER DECODER STAKE (DECRSTAKE-10)
 DRIVEN INTO VALVE BOX BASE
 GRAVEL/SOIL. SECURE DECODER TO
 STAKE WITH BOTTOM OF DECODER
 FACING UP AND FASTENED TO STAKE
 USING TWO (2) PROVIDED ZIP TIES
- 3 INCOMING HUNTER JACKETED IDWIRE PATH PER PLAN WITH 36-IN [1m] LOOP OF SLACK WIRE, MEASURED FROM TOP OF VALVE BOX, NEATLY COILED INSIDE VALVE BOX
- 4 RED AND BLUE WIRES FROM DECODER TO IDWIRE SPLICES
- 5 BARE COPPER GROUND WIRE FROM DECODER TO GROUNDING DEVICE, ROD OR PLATE, SEE PLAN FOR GROUND LOCATIONS AND TYPE
- 6 STATION WIRE FROM DECODER OUTPUT TO VALVE SOLENOID WIRES
- 7 IDWIRE PATH SPLICE WITH DECODER
 OUTPUT USING 3M DBRY/6 INSTALLED PER
 MANUFACTURERS INSTRUCTIONS
- 8 VALVE SOLENOID WIRE SPLICE WITH DECODER OUTPUT USING 3M DBRY/6 INSTALLED PER MANUFACTURERS INSTRUCTIONS
- 9 VALVE SOLENOID WIRE WITH 18-IN [45cm] COIL TO SPLICE WITH DECODER OUTPUT
- 10 STATION WIRE FROM DECODER OUTPUT TO VALVE SOLENOID WIRES IN A SEPARATE VALVE BOX
- (11) WIRE CONNECTION FROM DECODER OUTPUT TO SEPARATE SOLENOID SHALL BE 14 AWG, PE INSULATED

- OUTGOING HUNTER JACKETED IDWIRE PATH PER PLAN WITH 36-IN [1m] LOOP OF SLACK WIRE, MEASURED FROM TOP OF VALVE BOX, NEATLY COILED INSIDE VALVE BOX
- 13 IRRIGATION VALVE, VALVE BOX, AND RELATED ITEMS PER IRRIGATION PLAN AND DETAILS

DETAIL NOTES:

- A. CONTRACTOR SHALL INDICATE ASSOCIATED VALVE NUMBER(S) ON MANUFACTURER PROVIDED LABEL ON DECODER WITH PERMANENT MARKER
- B. CONTRACTOR SHALL LOCATE SPLICE WITH IDWIRE AND DECODER OUTPUT IN THE MIDDLE OF 36-IN[1m] ID WIRE LOOP IN VALVE BOX
- C. CONTRACTOR SHALL CONNECT BARE COPPER GROUND WIRE FROM DECODER OUTPUT TO GROUND DEVICE WIRE USING ERICO PG11L KIT PER MANUFACTURES INSTRUCTIONS
- D. 14 AWG PE WIRE DISTANCE FROM OUTPUT SOLENOID TO SOLENOID IN SEPARATE VALVE BOX SHALL NOT EXCEED 150-FT [45m]
- E. EARTH GROUND TO BE A MINIMUM OF 8 FEET AWAY FROM DECODER AND AT A RIGHT ANGLE TO THE TWO WIRE PATH
- F. GROUNDING SHALL BE INSTALLED AT A MINIMUM OF EVERY 12^{TH} DECODER OR 1,000', WHICHEVER OCCURS FIRST





ICD-200 DECODER